



sartorius

arium[®] pro Ultrapure Water Systems

Application-orientated and flexible to meet the highest demands



arium[®] pro ultrapure water systems

■ Modular – selection of five systems specially for your application

■ Flexible – perfect fit to every laboratory environment

■ Easy to use – glass touch display with intuitive menu navigation

Description

The arium[®] pro series offers a flexible and modular system which, compared to conventional devices, demonstrates excellent cost efficiency.

All systems meet and exceed the ASTM Type 1 water quality standards and ensure the best reproducible results in their class. Up to 120 l of consistently high-quality ultrapure water with a conductivity of 0.055 $\mu\text{S}/\text{cm}$ ($\approx 18.2 \text{ M}\Omega \times \text{cm}$) can be dispensed each hour. The ultrapure water is virtually microorganism-free when a Sartopore[®] 2 150 end filter is used.

The patented Sartorius technology, SD card slot, long service life, and low maintenance requirement distinguish the arium[®] pro systems as easy-to-use, efficient and reliable Type 1 ultrapure water systems.

Applications

- HPLC, GC-MS
- Biological applications e.g. Cell Culture
- Standard Applications such as preparation of solutions, media and reagents

Modular

The selection of five systems contains module components specially tailored to your application. arium[®] pro supplies the desired level of ultrapure water quality for any critical and standard applications.


Flexible

The space-saving installation of the device on, under, or above your workstation integrates it perfectly into any laboratory. The display-dispense unit can be flexible positioned.

Display with touch-function

Even the opening of the dispensing valve can be controlled by the unique display. Simply navigate intuitively in the easy-to-use and clear menu navigation by lightly touching the display – even with gloves.

Technical Specifications



Flow chart arium® pro VF TOC

Technical Specifications

Dimensions:	35.0 × 49.2 × 45,1 cm Width × height × depth
Empty weight	17 – 19 kg, depending on system type
Operating weight	27 – 29 kg, depending on system type
Power supply	100 – 240 VAC (± 10%); 50 – 60 Hz, 130 VA (max.)
Operating temperature	2 °C – 35 °C at max. 80% rel. humidity
Storage temperature	5 °C – 45 °C at max. 80% rel. humidity
Data output	SD card slot ² , RS-232 interface

Feed Water Quality


Water purified by reverse osmosis,
distillation or Deionisation.¹

Input pressure	0 – 6.9 bar, recomb. > 2 bar
Temperature	2 – 30 °C
Specific conductivity	< 100 µS/cm compensated to 25 °C
TOC content	< 50 ppb
Turbidity	< 1 NTU
pH value	4 – 10

¹ To operate arium® pro with non-treated drinking water the Universal Kit could be used in most cases. In order to verify the specifications of your feed water, please contact the Sartorius Application Support.

² Not applicable for arium® pro

arium® pro DI and arium® pro



Standard Applications

- AAS, ICP-MS
- Ion chromatography
- Preparation of reagents
- Photometry

Description

The arium® pro DI is a highly efficient water treatment system and the ultrapure water quality exceeds the ASTM Type 1 quality standard.

Water is purified by a three-stage process. In the first two stages, both organic and inorganic components are removed reliably from the feed water by using the specially developed Elemental Kit cartridges. The third purification stage is performed using a Sartopore® final filter, which is connected directly at the point of use and removes particles and bacteria.

The arium® pro offers an even more affordable alternative. Reduced to the essential features, it produces ultrapure water – without any compromises.

Product water quality

Water purification methods	Adsorption by spherical activated carbon, deionization, optional final particle and sterile filtration
Type of water	ASTM Type 1 ultrapure water
Output	120 l/h
Water dispensing flow rate	0.1 – 2 l/ min. adjustable
Volume controlled output ⁴	2 l/min in 100 ml, 1 l or 5 l steps, depending on total dispense between 0,1 l and 60 l
Volume accuracy ⁵	3 % in the range of 0.25 l and 60 l
Conductivity ¹	0.055 µS/cm compensated to 25 °C
Resistivity ¹	18.2 MΩ × cm compensated to 25 °C
TOC content ³	≤ 5 ppb
Microorganism content ²	< 1 CFU/1,000 ml
Particle content ²	< 1/ml

arium® pro DI and pro systems for producing ASTM Type 1 ultrapure water

Equipment supplied: 1 arium® pro and Connection Set

Order Number	Description
H2Opro-DI-T	arium® pro DI Bench-Top system in a compact design for every laboratory workstation
H2Opro-DI-B	arium® pro DI space saving Wall-Mounted system with integrated wall bracket
H2Opro-DI-D	arium® pro DI Built-In system, space saving in every laboratory cabinet, with integrated wall bracket for display-dispense unit
H2Obasic-T	arium® pro Bench-Top system in a compact design for every laboratory workstation
H2Obasic-B	arium® pro space saving Wall-Mounted system with integrated wall bracket

¹ Measured value output is adjustable to 25 °C compensated or non-compensated

² When using a Sartopore® 2 150 end filter

³ Feed water < 50 ppb TOC content

⁴ At 2 bar pressure, depending on the connected accessories and end filter

⁵ Under constant operating conditions

arium® pro UV



Chemical-Analytical Applications

- HPLC
- GC-MS, AAS, ICP-MS
- Ion chromatography
- TOC-Analysis
- Photometry

Description

Like arium® pro DI, the arium® pro UV features three-stage purification technology. Yet it additionally uses photooxidation to remove organic components.

With two different wavelengths, the horizontally positioned UV lamp (185 | 254 nm) prevents microbial growth and reliably reduces organic compounds to a TOC value of < 2 ppb.

Perfectly matched to support photooxidation technology, the Analytical Kit cartridges optimize water purification and specifically remove inorganic and organic substances.

The current TOC values are continuously measured by the optionally integrated TOC monitor and shown on the display.

Product water quality

Water purification methods	Adsorption by spherical activated carbon, de-ionization, UV-irradiation, optional final particle and sterile filtration
Type of water	ASTM Type 1 ultrapure water
Output	120 l/h
Water dispensing flow rate	0.1 – 2 l/ min. adjustable
Volume controlled output ⁴	2 l/min in 100 ml, 1 l or 5 l steps, depending on total dispense between 0,1 l and 60 l
Volume accuracy ⁵	3 % in the range of 0.25 l and 60 l
Conductivity ¹	0.055 µS/cm compensated to 25 °C
Resistivity ¹	18.2 MΩ × cm compensated to 25 °C
TOC content ³	< 2 ppb
Microorganism content ²	< 1 CFU/1.000 ml
Particle content ²	< 1/ml

arium® pro UV systems for producing ASTM Type 1 ultrapure water

Equipment supplied: 1 arium® pro with UV Lamp (185 | 254 nm) and Connection Set

Order Number	Description
H20pro-UV-T	arium® pro UV Bench-Top system in a compact design for every laboratory workstation, incl. UV Lamp
H20pro-UV-B	arium® pro UV space saving Wall-Mounted system with integrated wall bracket and incl. UV Lamp
H20pro-UV-D	arium® pro UV Built-In system, space saving in every laboratory cabinet, with integrated wall bracket for display-dispense unit and incl. UV Lamp
H20pro-UV-T-TOC	arium® pro UV Bench-Top system in a compact design for every laboratory workstation, incl. UV Lamp and TOC Monitor
H20pro-UV-B-TOC	arium® pro UV space saving Wall-Mounted system with integrated wall bracket, incl. UV Lamp and TOC Monitor
H20pro-UV-D-TOC	arium® pro UV Built-In system, space saving in every laboratory cabinet, with integrated wall bracket for display-dispense unit, incl. UV Lamp and TOC Monitor

¹ Measured value output is adjustable to 25 °C compensated or non-compensated

² When using a Sartopore® 2 150 end filter

³ Feed water < 50 ppb TOC content

⁴ At 2 bar pressure, depending on the connected accessories and end filter

⁵ Under constant operating conditions

arium® pro UF



Biological Applications

- AAS, ICP-MS
- Ion chromatography
- Electrophoresis
- Endotoxin-Analysis
- Immunocytochemistry
- Nutrient media for cell culture
- Production of monoclonal antibodies
- Photometry

Description

In addition to featuring highly efficient purification components, such as active carbon and highly effective ion exchange resin, the arium® pro UF ultrapure water system additionally includes a hollow-fiber ultrafilter. This ultrafilter uses crossflow technology to reliably remove endotoxins, microorganisms and particles, as well as DNases and RNases from ultrapure water. As a result, this filter makes the arium® pro UF ideal for use in cell cultivation, electrophoresis, and many other related applications.

Supported by the top-down flow technology incorporated into Biological Kit cartridges, the system produces ASTM Type 1 ultrapure water of the highest quality.

Product water quality

Water purification methods	Adsorption by spherical activated carbon, deionization, ultrafiltration, optional final particle and sterile filtration
Type of water	ASTM Type 1 ultrapure water
Output	120 l/h
Water dispensing flow rate	0.1 – 1.7 l/min. adjustable
Volume controlled output ⁴	1.7 l/min in 100 ml, 1 l or 5 l steps, depending on total dispense between 0,1 l and 60 l
Volume accuracy ⁵	3 % in the range of 0.25 l and 60 l
Conductivity ¹	0.055 µS/cm compensated to 25 °C
Resistivity ¹	18.2 MΩ × cm compensated to 25 °C
TOC content ³	≤ 5 ppb
Microorganism content ²	< 1 CFU/1,000 ml
Particle content ²	< 1/ml
Endotoxin	< 0.001 EU/ml
RNase content	< 0.004 ng/ml
DNase content	< 0.024 pg/µl

arium® pro UF systems for producing ASTM Type 1 ultrapure water

Equipment supplied: 1 arium® pro with Ultrafilter and Connection Set

Order Number	Description
H2Opro-UF-T	arium® pro UF Bench-Top system in a compact design for every laboratory workstation, incl. Ultrafilter
H2Opro-UF-B	arium® pro UF space saving Wall-Mounted system with integrated wall bracket and incl. Ultrafilter
H2Opro-UF-D	arium® pro UF Built-In system, space saving in every laboratory cabinet, with integrated wall bracket for display-dispense unit and incl. Ultrafilter

¹ Measured value output is adjustable to 25 °C compensated or non-compensated


² When using a Sartopore® 2 150 end filter

³ Feed water < 50 ppb TOC content

⁴ At 2 bar pressure, depending on the connected accessories and end filter

⁵ Under constant operating conditions

arium® pro VF



Chemical-Analytical, Biological, and Standard Applications

- HPLC, GC-MS, AAS, ICP-MS, IC
- PCR
- Electrophoresis
- Endotoxin-Analysis
- Immunocytochemistry
- Nutrient media for cell culture
- Production of monoclonal antibodies
- Photometry

Description

This high-end unit delivers ultrapure water of the highest quality and combines the purification technology of the arium® pro UF and pro UV units all into one system.

The integrated, horizontal UV lamp (185 | 254 nm), together with a hollow-fiber ultrafilter, not only prevents microbial growth and reduces the TOC content to a minimum of < 2 ppb, but also additionally removes endotoxins, microorganisms and particles, as well as DNases and RNases. As a result, it is the perfect solution for a large number of critical applications in the laboratory.

Current TOC values are continuously measured with the highest accuracy by the optionally integrated TOC monitor and displayed on the screen.

Product water quality

Water purification methods	Adsorption by spherical activated carbon, de-ionization, ultrafiltration, UV irradiation, optional final particle and sterile filtration
Type of water	ASTM Type 1 ultrapure water
Output	120 l/h
Water dispensing flow rate	0.1 – 1.7 l/min. adjustable
Volume controlled output ⁴	1.7 l/min in 100 ml, 1 l or 5 l steps, depending on total dispense between 0,1 l and 60 l
Volume accuracy ⁵	3 % in the range of 0.25 l and 60 l
Conductivity ¹	0.055 µS/cm compensated to 25 °C
Resistivity ¹	18.2 MΩ × cm compensated to 25 °C
TOC content ³	< 2 ppb
Microorganism content ²	< 1 CFU/1,000 ml
Particle content ²	< 1/ml
Endotoxin	< 0.001 EU/ml
RNase content	< 0.004 ng/ml
DNase content	< 0.024 pg/µl

arium® pro VF systems for producing ASTM Type 1 ultrapure water

Equipment supplied: 1 arium® pro with UV lamp (185 | 254 nm), Ultrafilter and Connection Set

Order Number	Description
H2Opro-VF-T	arium® pro VF Bench-Top system in a compact design for every laboratory workstation, incl. UV Lamp and Ultrafilter
H2Opro-VF-B	arium® pro VF space saving Wall-Mounted system with integrated wall bracket, incl. UV Lamp and Ultrafilter
H2Opro-VF-D	arium® pro VF Built-In system, space saving in every laboratory cabinet, with integrated wall bracket for display-dispense unit, incl. UV Lamp and Ultrafilter
H2Opro-VF-T-TOC	arium® pro VF Bench-Top system in a compact design for every laboratory workstation, incl. UV Lamp, Ultrafilter and TOC Monitor
H2Opro-VF-B-TOC	arium® pro VF space saving Wall-Mounted system with integrated wall bracket, incl. UV Lamp, Ultrafilter and TOC Monitor
H2Opro-VF-D-TOC	arium® pro VF Built-In system, space saving in every laboratory cabinet, with integrated wall bracket for display-dispense unit, incl. UV Lamp, Ultrafilter and TOC Monitor

¹ Measured value output is adjustable to 25 °C compensated or non-compensated

² When using a Sartopore® 2 150 end filter

³ Feed water < 50 ppb TOC content

⁴ At 2 bar pressure, depending on the connected accessories and end filter

⁵ Under constant operating conditions

arium® dispense gun

Ergonomic water dispensing with an up to 3.7 m working radius



- Extended working area of 3.7 m
- Available with a height-adjustable stand or wall mounting bracket
- Ergonomic design
- Easy-to-use
- Sterile-grade filters connectable

Description

The arium® dispense gun is an ergonomic, easy-to-handle dispensing unit that is excellently suited for dispensing ultrapure water.

Depending on the working environment, you can save space by mounting the dispense gun on the wall or on a stand that is height-adjustable up to 70 cm.

Moreover, the stand lets you work in a relaxed position and allows optimal adjustment to different sized sampling vessels. Thanks to the extended tubing, the work area has an expanded radius to up to 2.5 m away from the arium® unit and another 1.2 m to the stand.

It is also easy to install a sterile-grade filter Sartopore® 2 150 capsule with 0.2 µm pore size for guaranteed sterile and particulate-free water dispensing.

Technical Specifications

Materials

Stand	Aluminum (gray anodized)
Dispense gun	Plastic, white finish
Tubing	PVDF

Dimensions without tubing [W × H × D]

Dispense gun with stand	18.5 × 59.5 × 51.0 cm
Dispense gun with wall mounting bracket	9.0 × 10.0 × 28.5 cm

Weight without tubing

Dispense gun with stand	5.60 kg
Dispense gun with wall mounting bracket	0.46 kg

Intended Use

arium® comfort I and comfort II
arium® pro DI, pro UF, pro UV and pro VF
arium® 611

Order Number	Description
--------------	-------------

H2Opro-AMDG1	arium® Dispense Gun inclusive height-adjustable Stand, qty. 1 unit
H2Opro-AMDG2	arium® Dispense Gun inclusive wall mounting kit, qty. 1 unit

Multifunction stand

All menu functions directly at the dispensing port



- Visual quality control directly at the water-dispensing port
- System control directly at the work place
- Optimal customization to the variously sized sampling vessels
- Water dispensing using the slider
- Radius extends up to 2.5 meters
- Sterile-grade filters connectable

Description

The multifunction stand not only extends the working radius up to 2.5 meter from arium® ultrapure water system, but at the same time enables the control of the unit and monitors the ultrapure water quality directly at the dispensing port. The bracket is integrated in the stand to which the flexible arium® display is mounted. This creates a combination that provides complete access to the menu with all its functions along with the practical convenience of a dispense gun.

Moreover, the over 70 cm height-adjustable stand lets you work in a relaxed position one-handedly and allows optimal adjustment to different sized sampling vessels.

It is also easy to install a sterile-grade filter Sartopore® 2 150 capsule with 0.2 µm pore size for guaranteed sterile and particulate-free water dispensing.

Technical Specifications

Materials

Stand	Aluminum (gray anodized)
-------	-----------------------------

Dimensions [W × H × D]	22.0 × 59.5 × 25.5 cm
-----------------------------------	-----------------------

Weight	5.60 kg
---------------	---------

Intended use for arium® built-in units (D-version):

arium® comfort I and comfort II
arium® pro DI, pro UF, pro UV and pro VF

Order Number

Description

H2Opro-ADM1

arium® multifunction stand, height-adjustable, including display-mounting kit, for arium® Built-In systems. qty. 1 unit

Water Guard

Early Detection of Leakages Protects the Laboratory



- Highly sensitive optical sensor
- Audiovisual alarm signals
- Automatic water-stop in case of leakage
- High-grade material, non-corrosive
- Easy to install
- Integrated wall mounting bracket for magnetic valve

Description

Only the early detection of water leakages provides optimal protection against water damage in the laboratory. Leakages are registered by the highly sensitive optical sensor.

In contrast to conventional sensors, this sensor functions independently of conductivity measurement values as these are so low in the ultrapure water area that the activation of the guard is not guaranteed. Once a leakage is detected the water guard automatically locks the feed water inlet line. An acoustic warning is triggered immediately and the system status can be constantly controlled using the integrated LED display. The sensitive, optical sensors and high-grade materials mean that the arium® water guard is perfect for all ultrapure and pure water systems.

Technical Specifications

Sensor Dimensions

Diameter	5 cm
Height	2.5 cm
Cable length	2 m

Tubing Connections

Input	3/8" plug connection
Output	3/8" plug connection
Power supply	100 – 240 VAC 50 – 60 Hz

Intended Use

System type:
arium® comfort I and comfort II
arium® pro, pro DI, pro UF, pro UV and pro VF
arium® advance RO and EDI
arium® 611, 612 and 613

Order Number

Description

610AWG1	arium® water guard, qty. 1 unit
---------	---------------------------------

Foot switch

More convenient pure water dispensing



- Water dispensing at a press of the foot
- Ideally suited for clean rooms, minimize the risk of contamination
- Comfortable fatigue-free switching thanks to low height

Description

Easy-to-connect foot switch for starting and stopping the water dispensing. The rugged foot switch frees up both hands for other tasks, such as changing vessels, and minimizes the risk of contamination in the cleanroom.

Technical Specifications

Material	Nylon, glass fiber reinforced
Dimensions [W × H × D]	14.0 × 4.5 (max.) × 10.6 cm
Cable length	2 m
Power supply	100 – 240 VAC 50 – 60 Hz
Connector	Phoenix plug, 2-pin

Intended Use

arium® comfort I and comfort II
arium® pro DI, pro UF, pro UV and pro VF

Order Number

Description

H20-AFS1	arium® foot switch, qty. 1 unit
----------	---------------------------------

Level Sensor

Practically separate tank filling



- Flexible water transport to anywhere
- Can fill every tank system

Description

The level sensor makes it easy to connect an external water storage tank and subsequently fill a tank with ultrapure water.

Technical Specifications

Length Level Sensor	88 mm
Diameter Connection	2.03 cm (max.)
Bore	1.65 cm
Cable length	3 m

Intended Use

arium® pro DI, pro UF, pro UV and pro VF

Order Number

Description

H20-ALS1	arium® Level Sensor, qty. 1 unit
----------	----------------------------------

Printer

GMP Data logging make easy



- Acquisition and documentation of current measurement data

Description

Current measurements can be output to the printer via an RS-232 interface to support qualification and documentation tasks.

Technical Specifications

Dimensions [W×H×D]	21.5×7.8×15.5 cm
--------------------	------------------

Intended Use

arium® comfort I and comfort II
arium® pro DI, pro UF, pro UV and pro VF
arium® 611, 612 and 613

Order Number

Description

611APR1	arium® printer, qty. 1 unit
---------	-----------------------------

arium® pro Cartridge Sets

Pre-treatment and Post-treatment Cartridge Utilizing Top-Down Technology



- High performance capacity thanks to efficient ion-exchange resin
- Fast and effective absorption of impurities through high-grade activated carbon
- Optimized crossflow behavior, prevents separation of the resin mixed-bed
- Patented connection process simplifies the replacement of consumables

Description

The cartridge sets are optimized for the removal of both organic and inorganic constituents. Every set has been designed specifically to match the unit and delivers ultrapure water that even exceeds the ASTM type 1 quality standard. Such consistently high water quality is a guarantee for optimal reproducibility of your results.

Optimized cartridge materials such as highly effective activated carbon coupled with a highly efficient ion-exchange resin deliver long lasting performance and thereby ensure long maintenance intervals.

The Top-Down-Flow technology produces ideal purification kinetics and prevents any mixing of cleaning media. The cartridge was designed with the applicable standards for flow rate in the cross section and contact time with the medium in mind.

Technical Specifications

Materials

Housing	high-grade polypropylene
Mounting screws	stainless steel
Cleaning media	spherical catalytic effective activated carbon ultrapure, mixed bed ion exchange resin
Feed water requirements	see "Technical specifications" page 2

Exchange capacity at 18.2 MΩ × cm ultrapure water related to CaCO₃

	[Grain]	[Equivalent]
Analytical Kit	965	1.25
Biological Kit	1,141	1.48
Elemental Kit	1,268	1.64
Universal Kit	965	1.25

Intended Use

H2O-A-PACK	arium® pro VF and pro UV
H2O-B-PACK	arium® pro UF
H2O-E-PACK	arium® pro and pro DI
H2O-U-PACK*	arium® pro, pro DI, pro UF, pro UV and pro VF

Order Number	Description
H2O-A-PACK	Analytical Kit, arium® pro Cartridge Set for biological, chemical-analytical and standard ultrapure water applications, qty. 1 unit
H2O-B-PACK	Biological Kit, arium® pro Cartridge Set for biological ultrapure water applications, qty. 1 unit
H2O-E-PACK	Elemental Kit, arium® pro Cartridge Set for standard ultrapure water applications, qty. 1 unit
H2O-U-PACK*	Universal Kit, arium® pro Cartridge Set for non-treated drinking water*, qty. 1 unit

*) To operate arium® pro with non-treated drinking water the Universal Kit could be used in most cases. In order to verify the specifications of your feed water, please contact the Sartorius Application Support.

Sterile-grade Sartopore® 2 150 Capsules

Sterile and Particle-free Water Dispensing



- Excellent service lifetime and flow rate
- Integrity tested
- Validated acc. to HIMA and ASTM F-838-05
- Meets WFI quality standards acc. to USP incl. USP plastic class VI test
- Manufacture acc. to DIN ISO 9001
- Easy to install
- Automatic venting
- Certified quality

Description

Sartopore® 2 150 is a sterile, ready-to-use membrane filter capsule for the most stringent requirements. Sartopore® 2 150 membrane filter capsules contain a hydrophilic, heterogeneous polyethersulfone double membrane. It enables an excellent service life and output. The capsule is attached, by a quick connector, at the final position and reliably removes all particles. A hydrophobic PTFE membrane at the farthest point „upstream“ allows for easy and clean venting of the capsule.

All pleated Sartopore® 2 membrane filter units are validated as sterile filters for biopharmaceutical applications according to the HIMA and ASTM F-838-05 guidelines (documentation available). During the manufacturing process, every unit is tested for integrity in order to meet the highest quality standards and safety regulations.

Technical Specifications

Materials

Membranes	Asym. Polyethersulfone
Filling bell	Polycarbonate
Other plastics	Polypropylene
Pore size	0.45 µm + 0.2 µm
Filtration area	0.015 m ²
Input and output	1/4" plug connection
Sterilization (max. 3 cycles)	Autoclaving at 134°C, 1 bar, 30 min.
Max. diffusion	1 ml/min @ 2.5 bar
Min. bubble point	3.2 bar

Intended Use

On Dispense Guns and Display-Dispense Unit for system type:

arium® comfort I and comfort II
arium® pro, pro DI, pro UF, pro UV and pro VF
arium® 611
arium® bagtank Dispense Gun
arium® Dispense Gun

Order Number	Description
--------------	-------------

5441307H4--CE--B	Sterile-grade Sartopore® 2 150 capsules, 0.2 µm pore size, qty. 5 units
------------------	---

UV Lamp (185 | 254 nm)

Ultrapure water, free of TOC



- Horizontal installation, optimize temperature gradient
- Effectively destroys organic compounds
- Reduces microbiological growth
- Easy replacement

Description

The horizontally arranged UV lamp delivers especially reliable results. Unlike vertical units, the temperature gradient is less pronounced and does not affect the activity of UV waves. The two different wavelengths reliably removes organic substances (TOC or total organic carbon), effectively preventing microbiological growth. Organic compounds oxidize at 185 nm, whereas 254 nm prevents microbiological growth.

Technical Specifications

Material	quartz glass
----------	--------------

TOC-content product water*	< 2 ppb
----------------------------	---------

Intended Use

arium® comfort I and comfort II
arium® pro UV and pro VF
arium® 611 UV and VF

* Feed water < 50 ppb TOC content

Order Number

Description

611CEL1	arium® UV Lamp (185 254 nm), qty. 1 unit
---------	--

Ultrafilter

Ultrapure water, free of endotoxines, DNases and RNases



- High flow rates
- Integrity tested
- Long service lives
- Certified quality

Description

The hollow-fiber ultrafilter utilizes crossflow technology to reliably remove bacterial endotoxins, microorganisms and particulates, as well as DNases and RNases from the ultrapure water.

The filters have been developed and manufactured according to DIN EN ISO 9001 | DIN EN 46 001 certified quality assurance system that meets the requirements of the FDA's Quality System Regulation 21 CFR Part 820. During the manufacturing process, every unit is integrity tested to meet the highest quality standards and safety regulations.

Technical Specifications

Materials

Membrane	polysulfone
Composites	polyurethane (PUR)
Housings, Caps	polycarbonate (PC)
Plugs	polypropylene (PP)
Effective membrane area	2.1 m ²
Max. pressure	3 bar at room temperature
Sanitization	200 ppm sodium hypochlorite, 45 min, max. 1x/week

Filling volume

Lumen	152 ml
Filtrate area	306 ml

Retention of bacteria and endotoxins

<i>Brev. diminuta</i>	LRV 7-10
<i>E. coli</i> O55:B5 Endotoxin	LRV > 3.5
Natural Endotoxins	LRV > 3.0

Dimensions fibers

Inner diameter	215 µm
Wall thickness	50 µm
Molecular Weight Cut Off (MWCO)	5.000 (= 5 kD)

Intended Use

arium[®] pro VF and pro UF
arium[®] 611 VF and UF

Order Number	Description
611CDU5	arium [®] Ultrafilter, qty. 1 unit

Biofilm Cleaning Syringes

Effective removal of microorganisms for a long lifetime



- Highly effective against biofilms (consisting of bacteria, fungi etc.)
- Free of organic components (TOC)
- Tenside-free
- Gentle on the environment and materials

Description

With this cleaning agent, the regular removal of the biofilms that develop during water purification is easy and effective. Biofilms especially develop on the concentrate side of ultrafilters. If biofilms are removed regularly, higher flow rates and longer lifetimes result.

This highly effective solution is prefilled in 50 ml syringes and directly ready for use. It does not form any trihalogenmethane, chloramines, hydrochloric acid or chlorate. The cleaning agent is non-caustic and degrades into NaCl and water.

Technical Specifications

Ingredients

- Sodium hypochlorite
- Hydrogen peroxide
- Sodium chloride
- Hypochlorous acid

Intended Use

arium® pro DI, pro UF, pro UV and pro VF
arium® 611

Order Number

Description

611CDS2	arium® Biofilm Cleaning Syringes, prefilled in 50 ml syringes, qty. 2 units
611CDS6	arium® Biofilm Cleaning Syringes, prefilled in 50 ml syringes, qty. 6 units

Sales and Service Contacts

For further contacts, visit www.sartorius.com

Europe

Germany

Sartorius Weighing Technology GmbH
Weender Landstrasse 94-108
37075 Goettingen
Phone +49.551.308.0
Fax +49.551.308.3289
info.mechatronics@sartorius.com
www.sartorius.com

Austria

Sartorius Mechatronics Austria GmbH
Franzosengraben 12
1030 Wien
Phone +43.1.7965760.0
Fax +43.1.7965760.24
info.austria@sartorius.com

France Et Suisse Romande

Sartorius Mechatronics
France SAS
3 Avenue du Canada
91940 Les Ulis
Phone +33 (0) 1 69 19 21 00
Fax +33 (0) 1 69 20 09 22
service.client@sartorius.com

Belgium

Sartorius Mechatronics
Belgium N.V.
Leuvensesteenweg, 248/B
1800 Vilvoorde
Phone +32.2.756.06.71
Fax +32.2.253.45.95
info.belgium@sartorius.com

Hungary

Sartorius Mechatronics Hungária Kft.
Kagyló u. 5.
2092 Budakeszi
Phone +3623.457.227, 457.228, 457.148
Fax +3623.457.147
mechatronika@sartorius.hu

Ireland

Sartorius Mechatronics Ireland Limited
Unit 41, The Business Centre
Stadium Business Park
Ballycoolin Road
Dublin 11
Phone +353-(0)1-8089050
Fax +353-(0)1-8089388
info.ireland@sartorius.com

Italy

Sartorius Mechatronics Italy S.r.l.
Uffici di Milano
Viale A. Casati, 4
20053 Muggiò (Milan)
Phone +39.039.46591
Fax +39.039.465988
info@sartorius.it

Netherlands

Sartorius Mechatronics
Netherlands B.V.
Edisonbaan 24
3439 MN Nieuwegein
Phone +31.30.6053001
Fax +31.30.6052917
weegtechnik.nl@sartorius.com

Poland

Sartorius Mechatronics
Poland Sp. z o.o.
ul. Wrzesinska 70
62-025 Kostrzyn
Phone +48.61.6473830
Fax +48.61.6473839
info.pl@sartorius.com

Spain

Sartorius Mechatronics Spain S.A.U.
Offices in Madrid:
c/ Isabel Colbrand, 10-12, of. 70
28050 Madrid
Phone +34.91.358.60.94
Fax +34.91.358.84.85
Sartorius Mechatronics Spain S.A.U.
Offices in Barcelona:
C/Marcus Porcius, 1 (Edificio BCIN)
Poligon Les Guixeres s/n
08915 - Badalona
Barcelona - Spain
Phone +34.902.123.367
Fax +34.91.358.96.23
spain.weighing@sartorius.com

Switzerland

Sartorius Mechatronics Switzerland AG
Ringstrasse 24a
8317 Tagelswangen (ZH)
Phone +41.44.746.50.00
Fax +41.44.746.50.50
mechatronics.switzerland@sartorius.com

U.K.

Sartorius Mechatronics UK Ltd.
Longmead Business Centre
Blenheim Road, Epsom
Surrey. KT19 9QQ
Phone +44.1372.737102
Fax +44.1372.729927
uk.customerservice@sartorius.com

America

Argentina

Sartorius Argentina S.A.
Int. A. Ávalos 4251
B1605ECS Munro
Buenos Aires
Phone +54.11.4721.0505
Fax +54.11.4762.2333
sartorius@sartorius.com.ar

Brazil

Sartorius do Brasil Ltda.
Av. D. Pedro I, 241
Vila Pires - Santo André
São Paulo
09110-001
Phone +55.11.4451.6226
Fax +55.11.4451.4369
sartorius@sartorius.com.br

Canada

Sartorius Mechatronics Canada
2179 Dunwin Drive #4
Mississauga, ON L5L 1X2
Phone +1.905.569.7977
Toll-Free +1.800.668.4234
Fax +1.905.569.7021
sales.canada@sartorius.com

Mexico

Sartorius de México S.A. de C.V.
Circuito Circunvalación Poniente
No. 149
53100, Satélite Estado de México,
México
Phone +5255.5562.1102
Fax +5255.5562.2942
sartorius@sartomex.com.mx

USA

Sartorius Mechatronics
Corporation
5 Orville Drive
Bohemia, NY 11716
Phone +1.631.254.4249
Toll-free +1.800.635.2906
Fax +1.631.254.4253
wt.sales@sartorius.com

Asia | Pacific

China

Sartorius Scientific
Instruments (Beijing) Co., Ltd.
Konggang Industrial Zone B
No. 33 Yu'an Road
101300 Beijing, Shunyi District
Phone +86.10.8042.6300
Fax +86.10.8042.6486
ssil@sartorius.com

Hong Kong

Sartorius Mechatronics
Hong Kong Ltd.
Unit 1012, Lu Plaza
2 Wing Yip Street
Kwung Tong
Kowloon, Hong Kong
Phone +852.2774.2678
Fax +852.2766.3526
enquiry.hongkong@sartorius.com

India

Sartorius Mechatronics India Pvt Ltd.
69/2 Et 69/3, Jakkasandra,
Kunigal Road, Nelamangala Tq
Bangalore-562 123
Phone +91.80.4350.5250/51/52
mechatronics-india@sartorius.com

Japan

Sartorius Mechatronics Japan K.K.
8-11, Kita-Shinagawa 1-chome
Shinagawa-ku
Tokyo 140-0001
Phone +81.3.3740.5408
Fax +81.3.3740.5406
info@sartorius.co.jp

Philippines

Sartorius Mechatronics Philippines,
Incorporated
Unit 20-A The World Centre Building
330 Senator Gil Puyat Avenue Makati
City Philippines 1209
Phone +632.8640929
Fax +632.8640932
enquiry.philippines@sartorius.com

Singapore

Sartorius Mechatronics
Singapore Pte. Ltd.
1 Science Park Road #05-08A
The Capricorn
Singapore Science Park II
Singapore 117528
Phone +65.6872.3966
Fax +65.6778.2494
enquiry.singapore@sartorius.com

South Korea

Sartorius Mechatronics
Korea Ltd.
Yangjae B/D 4, 5F
209-3, Yangjae-Dong, Seocho-Ku
137-893 Seoul, Korea
Phone +82.2.575.6945
Fax +82.2.575.6949
enquiry.korea@sartorius.com

Thailand

Sartorius Mechatronics
Thailand Co. Ltd.
No. 129 Rama IX Road
Huaykwang
Bangkok 10310
Phone +66 2643.8361
Fax +66 2643.8367
enquiry.thailand@sartorius.com

Australia

Sartorius Mechatronics
Australia Pty Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175
Phone +61.3.8762.1800
Fax +61.3.8762.1828
Info.Australia@Sartorius-Stedim.com



◀ www.sartorius.com